

# Chemistry 254

## Experiment 4

### Solubility of a solid

Adam Menne  
Stellenbosch University

Last edited on 2nd September 2022

#### **Abstract**

In this practical the temperature dependence on the solubility of benzoic acid is investigated. This is related to the enthalpy of solution by use of the Van 't Hoff equation

## **Contents**

|                                    |          |
|------------------------------------|----------|
| <b>1 Results</b>                   | <b>2</b> |
| <b>2 Discussion</b>                | <b>2</b> |
| <b>Appendix A Additional tasks</b> | <b>3</b> |
| <b>Appendix B Flow Diagram</b>     | <b>4</b> |
| <b>Appendix C MSDS</b>             | <b>5</b> |

## 1 Results

At 20 °C the concentration of benzoic acid was 0.02173 *M*, while at 30 °C it was 0.02929 *M*

The temperature independent  $\Delta H^0$  was calculated to be 22.064 *kJ*

A static export of the notebook containing all analysis and figures is available at [https://adammenne.github.io/chemistry\\_254/practical\\_4/notebook.html](https://adammenne.github.io/chemistry_254/practical_4/notebook.html).  
With full source code available at [https://github.com/AdamMenne/chemistry\\_254/tree/master/practical\\_4](https://github.com/AdamMenne/chemistry_254/tree/master/practical_4)

## 2 Discussion

The titrations carried out were consistent with both sets having standard deviations below 0.00014.

## Appendix A Additional tasks

1. 0.02529
2. This rearrangement of the Van 't Hoff equation is possible due to the fact that  $\frac{d\frac{1}{T}}{dT} = -\frac{1}{T^2}$
3. The Van 't Hoff plot for a reaction that has a temperature dependant standard enthalpy will not be linear. With an adequate sample size over an appropriate range of temperatures, a polynomial fit will allow the gradient at any temperature to be found, and thus the standard enthalpy of solution to be calculated.

## Appendix B    Flow Diagram

1. Heat a 25  $\text{cm}^3$  pipette with a hair dryer, attach cotton filter, and fill past the mark
2. Remove filter, wipe off crystals, and evacuate the pipette to the mark
3. Empty into Erlenmeyer flask, washing with ethanol, add 2 to 3 drops phenolphthalein
4. Titrate with 0.1  $M$   $\text{NaOH}$
5. Clean pipette with ethanol followed by acetone, and dry
6. Repeat 3 times for both temperatures

## Appendix C   MSDS

### Benzoic acid

- Corrosive, health hazard
  - causes serious eye damage
  - may cause respiratory irritation
  - if in contact with skin or eyes wash for several minutes

### Sodium hydroxide

- Corrosive
  - may cause skin burns, eye damage
  - if in contact with skin or eyes wash for several minutes